



重组人 FGD3 蛋白产品说明书

一、产品基础信息

货号: REP06881

蛋白名称: Fgd1 Family, Member Faciogenital Dysplasia Fyve, Rhogef And Ph Domain Containing Zinc Finger Fyve Domain-Containing Protein Fyve, Rhogef And Ph Domain-Containing Protein 3 (FGD3)

别名: FGD1 family, member faciogenital dysplasia FYVE, RhoGEF and PH domain containing zinc finger FYVE domain-containing protein FYVE, RhoGEF and PH domain-containing protein 3

UniProt 编号: Q5JSP0

UniProt 查询链接: <https://www.uniprot.org/uniprotkb/Q5JSP0/entry>

种属来源: 人 (Homo sapiens)

表达系统: 大肠杆菌 (E.coli)

蛋白片段: 全长 (1-725 位氨基酸)

分子量: 86.3 kDa

蛋白标签: C 端 6×His 标签

二、氨基酸序列 (1-725aa)

MESGRGSSTPPGPAAALGMPDTGPGSSSLGKLQALPVGPRAHCGDPVSLAAAGDG
SPDIGPTGELSGSLKIPNRDSGIDSPSSSVAGENFPCEEGLAAGPSPTVLGAHAEMA
LDSQVPKVTPQEEADSDVGEEPDSENTPQKADKDALAQHSGPQKLLHIAQELLHT
EETVVKRLHLLDQVFCTRLTDAGIPPEVIMGIFSNISSIHRFHGQFLLPELKTRITEEW
DTNPRLGDILQKLAPFLKMYGEYVKNFDRVGLVSTWTQRSPLFKDVVHSIQKQEV
CGNLTQLHHMLEPVQRVPRYELLLKDYLRPLQDAPDRKDAERSLELISTAANHSN
AAIRKVEKMHKLLEVYEQLGGEEDIVNPANELIKEGQIQKLSAKNGTPQDRHLFLFN
SMILYCVPKRLRMGQKFSVREKMDISGLQVQDIVKPNTAHTFIITGRKRSLELQTRTE
EEKKEWIIQATIEKHKQNSSETFKAFFGAFSQDEDPSPDMPITSTSPVEPVVTE
GSSGAAGLEPRKLSKTRRDKEKQCKSCGETFNSITKRRHHCKLCGAVICGKCSE
FKAENSRQSRVCRDCFLTQPVAPESTEKTPTADPQPSLLCGPLRLSESGETWSEV
WAAIPMSDPQVLHLQGGSQDGRPLPRTIPLPSCKLSVPDPEERLDSGHVWKLQWAK
QSWYLSASSAELQQQWLETLSTAAHGDTAQDSPGALQLQVPMGAAAP



三、缓冲液体系

液体剂型：基于 Tris/PBS 缓冲体系，含 5%-50%甘油。可根据客户需求定制甘油含量，下单时可备注特殊要求。

冻干粉剂型：冻干前缓冲液为含 6%海藻糖的 Tris/PBS 缓冲体系。

四、储存条件

产品收到后可置于-20℃或-80℃环境保存。为保障蛋白活性、避免反复冻融导致蛋白变性失活，建议分装保存，按需取用。

五、产品简介

本产品为大肠杆菌表达系统制备的重组人 FGD3 蛋白全长（1-725 位氨基酸）。蛋白携带 C 端 6×His 标签，便于蛋白纯化、检测与鉴定。FGD3 (FGD1 family, member faciogenital dysplasia FYVE, RhoGEF and PH domain containing zinc finger FYVE domain-containing protein FYVE, RhoGEF and PH domain-containing protein 3) 相关蛋白产品，适用于体外功能实验、蛋白互作研究、抗体制备及各类生物医学科研场景。

六、注意事项

- 严禁反复冻融产品，最大程度保证蛋白稳定性与生物活性。
- 若对缓冲液组分有特殊要求，请在下单时备注说明。
- 本产品仅用于科学研究，不可用于临床诊断、治疗及商业生产等用途。