



重组人 SIRT7 蛋白产品说明书

一、产品基础信息

货号: REP08237

蛋白名称: Sirtuin Sirtuin Type Sir2-Like Protein Sir2-Related Protein Type Regulatory Protein Sir2 Homolog Nad-Dependent Deacetylase Sirtuin-Nad-Dependent Protein Deacylase Sirtuin-Nad-Dependent Protein Deacetylase Sirtuin-Silent Mating Type Information Regulation 2, S.Cerevisiae, Homolog 7 (SIRT7)

别名: sirtuin sirtuin type SIR2-like protein sir2-related protein type regulatory protein SIR2 homolog NAD-dependent deacetylase sirtuin-NAD-dependent protein deacetylase sirtuin-NAD-dependent protein deacety...

UniProt 编号: Q9NRC8

UniProt 查询链接: <https://www.uniprot.org/uniprotkb/Q9NRC8/entry>

种属来源: 人 (Homo sapiens)

表达系统: 大肠杆菌 (E.coli)

蛋白片段: 全长 (1-400 位氨基酸)

分子量: 52.3 kDa

蛋白标签: N 端 10×His 标签、C 端 Myc 标签

二、氨基酸序列 (1-400aa)

MAAGGLSRSERKAAERVRLREEQQRERLRQVSRILRKAAAERSAEEGRLLAESA
DLVTELQGRSRRREGLKRRQEEVCDDPEELRGKVRELASAVRNAKYLVVYTGAGIS
TAASIPDYRGPNGVWTLQKGRSVSAADLSEAEP TLTHMSITRLHEQKLVQHVV SQ
NCDGLHLRSGLPRTAISELHG NMYIEVCTSCVPNREYVRVFDV TERTALHRHQ TGR
TCHKCGTQLRDTIVHFGERGTLGQPLNWEAATEAASRADTILCLGSS LKVLKKYPRL
WCMTKPPSRRPKLYIVNLQWTPKDDWAALKLHGKCDDVMRLLMAELGLEIPAYS R
WQDPIFSLATPLRAGEEGSHSRKSLCRSREEAPPGDRGAPLSSAPILGGWFGRGC
TKRTRKKVT

三、缓冲液体系



液体剂型：基于 Tris/PBS 缓冲体系，含 5%-50%甘油。可根据客户需求定制甘油含量，下单时可备注特殊要求。

冻干粉剂型：冻干前缓冲液为含 6%海藻糖的 Tris/PBS 缓冲体系。

四、储存条件

产品收到后可置于-20℃或-80℃环境保存。为保障蛋白活性、避免反复冻融导致蛋白变性失活，建议分装保存，按需取用。

五、产品简介

本产品为大肠杆菌表达系统制备的重组人 SIRT7 蛋白全长（1-400 位氨基酸）。蛋白携带 N 端 10×His 标签、C 端 Myc 标签，便于蛋白纯化、检测与鉴定。SIRT7 (sirtuin sirtuin type SIR2-like protein sir2-related protein type regulatory protein SIR2 homolog NAD-dependent deacetylase sirtuin-NAD-dependent protein deacetylase sirtuin-NAD-dependent protein deacetylase sirtuin-silent mating type information regulation 2, *S.cerevisiae*, homolog 7) 相关蛋白产品，适用于体外功能实验、蛋白互作研究、抗体制备及各类生物医学科研场景。

六、注意事项

- 严禁反复冻融产品，最大程度保证蛋白稳定性与生物活性。
- 若对缓冲液组分有特殊要求，请在下单时备注说明。
- 本产品仅用于科学研究，不可用于临床诊断、治疗及商业生产等用途。